

Leading the World in  
G.I. Disposable Probes,  
Supplies & Accessories



# Reflux Monitoring

Comfortable...easy...  
smooth...flexible...

## Ambulatory Reflux


Ambulatory reflux probes from Diversatek Healthcare are the perfect choice for every application. Our pH and Impedance/pH probes are available in three sizes to meet a variety of esophageal lengths. Small in diameter, our probes are made from specially-formulated polyurethane\*, which softens in response to body temperature to optimize ease of insertion and enhance patient comfort. Not only that, our probes are 100% guaranteed...we stand behind our products!

### ComforTEC® Z/pH Probe

#### Internal Reference, Non-infused.

Single-use, non-sterile

Part No.	Description	Qty/Case
ZAN-BD-31	≥ 18 cm Esophageal length, 6 impedance, 2 pH channels 0 & 15 cm, 6.4 FR / 2.13 mm.	1
ZAN-BG-44	≥ 18 cm Esophageal length, 6 impedance, 2 pH channels -15 & 0 cm, 6.4 FR / 2.13 mm (Gastric).	1
ZAN-BS-01	≥ 18 cm Esophageal length, 6 impedance, 1 pH channel, 6.4 FR / 2.13 mm.	1
ZAN-BT-60	≥ 18 cm Esophageal length, 6 impedance, 3 pH channels -15 & 0 & 15 cm, 6.4 FR / 2.13 mm (Gastric).	1
ZPN-BS-46	15–18 cm Esophageal length, 6 impedance, 1 pH channel, 6.4 FR / 2.13 mm.	1
ZPN-BS-50	15–18 cm Esophageal length, 6 impedance, 1 pH channel, 6.4 FR / 2.13 mm. Sensor spacing at 1.5 cm.	1
ZIN-BS-51	< 15 cm Esophageal length, 6 impedance, 1 pH channel, 6.4 FR / 2.13 mm.	1





Diversatek Healthcare has a complete line of innovative products for Total Reflux Monitoring including Zvu<sup>®</sup> Advanced GI Diagnostic Software, ZepHr<sup>®</sup> Impedance/pH Recorder and ComforTEC<sup>®</sup> Z/pH Probes.



## Reflux Monitoring

### ZepHr<sup>®</sup> Impedance/pH Recorder & Accessories

#### ZepHr<sup>®</sup> Impedance/pH Recorder & Accessories

Part No.	Product	Description	Qty/Case
Z07-2000B-B	ZepHr <sup>®</sup> Z/pH Recorder	ZepHr <sup>®</sup> ambulatory impedance/pH reflux monitoring recorder. Includes carrying case, two SD cards, and two AA batteries. ZepHr <sup>®</sup> has a limited two year warranty.	1
Z07-9000	ZepHr <sup>®</sup> System Kit	Includes an SD Card Reader, pH Calibration Tubes, and Buffer Solutions.	1
Z07-4053	ZepHr <sup>®</sup> Battery Cover		1
Z07-4055	ZepHr <sup>®</sup> Probe Cover, Smoke		1
Z07R-8002	ZepHr <sup>®</sup> Thermoform Pouch		1
Z07-1250	SD Flash Memory Card		1
A04-2000-S	ZepHr <sup>®</sup> AirFlo™ Sphincter Locator	LES locator for pH or Z/pH reflux testing. Infused catheters are required for use.	1
D82-3103	pH 1 Buffer Solution (500 mL)	For use with pH reflux testing system and catheters.	1
D82-3151	pH 4 Buffer Solution (500 mL)	For use with Diversatek Healthcare Z/pH or pH reflux testing system and catheters.	1
D82-3152	pH 7 Buffer Solution (500 mL)	For use with Diversatek Healthcare Z/pH or pH reflux testing system and catheters.	1
D82-3181	pH 1 Buffer Solution (250 mL)	For use with pH reflux testing system and catheters.	1
D82-3182	pH 4 Buffer Solution (250 mL)	For use with Diversatek Healthcare Z/pH or pH reflux testing system and catheters.	1
D82-3183	pH 7 Buffer Solution (250 mL)	For use with Diversatek Healthcare Z/pH or pH reflux testing system and catheters.	1
L89-5700-4	USB Multicard Reader Kit for ZepHr <sup>®</sup>		1
L89-6020-40	pH Calibration Tubes (Set of 3)	Includes pH 4 and pH 7 buffer tubes plus a water tube.	1
L89-6021-40	pH Calibration Buffer Tube, 40 cm	Includes both pH 4 and pH 7 labels.	1
MOUSEPAD-Z1	Diversatek Healthcare Mousepad		1
NG-0050	Nasal Tube Fastener	Adult size for reflux testing system.	50




## ComfoTEC® Z/pH Probe

### Internal Reference, Infused.

Single-use, non-sterile. Use with ZepHr® AirFlo™ Sphincter Locator


Part No.	Description	Qty/Case
ZAI-BD-31	≥ 18 cm Esophageal length, 6 impedance, 2 pH channels 0 & 15 cm, 6.9 FR / 2.3 mm.	1
ZAI-BS-01	≥ 18 cm Esophageal length, 6 impedance, 1 pH channel, 6.9 FR / 2.3 mm.	1
ZPI-BS-46	15–18 cm Esophageal length, 6 impedance, 1 pH channel, 6.9 FR / 2.3 mm.	1
ZAI-BG-44	≥ 18 cm Esophageal length, 6 impedance, 2 pH channels -15 & 0 cm, 6.9 FR / 2.3 mm (Gastric).	1
ZPI-BG-52	15–18 cm Esophageal length, 6 impedance, 2 pH channels -12 & 0 cm, 6.9 FR / 2.3 mm. Sensor spacing 1.5 cm (Gastric).	1



### ComfoTEC® LPR Probe. Internal Reference, Infused.

Single-use, non-sterile


Part No.	Description	Qty/Case
ZAI-BL-54	26–28 cm Esophageal length, 6.9 FR / 2.3 mm. To assess acid and nonacid reflux (Total LPR).	1
ZAI-BL-55	23–25 cm Esophageal length, 6.9 FR / 2.3 mm. To assess acid and nonacid reflux (Total LPR).	1
ZAI-BL-56	20–22 cm Esophageal length, 6.9 FR / 2.3 mm. To assess acid and nonacid reflux (Total LPR).	1



### External Reference, Non-infused.

Single-use, non-sterile


Part No.	Description	Qty/Case
ZAN-BS-01E	≥ 18 cm Esophageal length, 6 impedance, 1 pH channel, 6.4 FR / 2.13 mm.	1
ZPN-BS-46E	15–18 cm Esophageal length, 6 impedance, 1 pH channel, 6.4 FR / 2.13 mm.	1
ZIN-BS-45E	< 15 cm Esophageal length, 6 impedance, 1 pH channel, 6.4 FR / 2.13 mm.	1
ZAN-BG-44E	≥ 18 cm Esophageal length, 6 impedance, 2 pH channels -15 & 0 cm, 6.4 FR / 2.13 mm. Special Order (Gastric).	1
ZPN-BG-07E	15–18 cm Esophageal length, 6 impedance, 2 pH channels -12 & 0 cm, 6.4 FR / 2.13 mm. Special Order (Gastric).	1




### External Reference, Infused.

Single-use, non-sterile

Part No.	Description	Qty/Case
ZAI-BS-01E	≥ 18 cm Esophageal length, 6 impedance, 1 pH channel, 6.9 FR / 2.3 mm.	1
ZAI-BG-44E	≥ 18 cm Esophageal length, 6 impedance, 2 pH channels -15 & 0 cm, 6.9 FR / 2.3 mm. Special Order (Gastric).	1






## ComforTEC® PLUS pH Probe

### Internal Reference, Non-infused.

Single-use, non-sterile


Part No.	Description	Qty/Case
PHN10	Dual channel 10 cm spacing, 6.4 FR / 2.13 mm.	1
PHN10-P	Dual channel 10 cm spacing, 5 FR / 1.67 mm.	1
PHN15	Dual channel 15 cm spacing, 6.4 FR / 2.13 mm.	1
PHN21	Dual channel 21 cm spacing, 6.4 FR / 2.13 mm.	1
PHN5	Dual channel 5 cm spacing, 6.4 FR / 2.13 mm.	1
PHNS	Single channel, 6.4 FR / 2.13 mm.	1
PHNS-P	Single channel, 5 FR / 1.67 mm.	1



### Internal Reference, Infused.

Single-use, non-sterile

Part No.	Description	Qty/Case
PHI10	Dual channel 10 cm spacing, 6.9 FR / 2.3 mm.	1
PHI15	Dual channel 15 cm spacing, 6.9 FR / 2.3 mm.	1
PHI5	Dual channel 5 cm spacing, 6.9 FR / 2.3 mm.	1
PHIS	Single channel, 6.9 FR / 2.3 mm.	1



# Innovations in Clinical Education

## Diversatek University Online

Our online training platform contains free content on esophageal and anorectal manometric studies, as well as impedance/pH reflux monitoring studies. Included are tutorials providing step-by-step guidance to develop skills in data acquisition, study review and report generation. Simply go to **DiversatekHealthcare.com** to request log-in information.

## Denver Training Center

Our Technical Research & Training Center specializes in clinical education and personalized training to make learning easy and convenient. Our comprehensive clinical instruction provides users with the knowledge and skills necessary to effectively acquire and analyze high resolution impedance manometry studies, impedance/pH reflux monitoring studies and anorectal manometry studies. Email us at **clinicaleducation@diversatekhc.com** or visit us online to learn more about our Denver course offerings.

## Webinars

Diversatek Healthcare is proud to present a series of live, interactive discussions on topics related to esophageal function testing, impedance/pH reflux monitoring studies and anorectal manometry. Each webinar includes a didactic session followed by an open discussion. All webinars are recorded and posted to the Diversatek U online portal for easy reference. Access **DiversatekHealthcare.com** for upcoming webinar announcements.

## The Diversatek Healthcare Review

Diversatek Healthcare Review, our monthly e-newsletter, features what's new at Diversatek University along with up-to-date product information. Every issue also includes our Clinical Insights, providing educational tips for Z/pH and HRiM analysis as well as answers to the most frequently asked questions.



# Personalized Clinical Support

## Onsite Training

Diversatek Healthcare Clinical Specialists deliver product support to suit your specific needs—on your schedule. Specialists are onsite at your facility to train and support you on your Diversatek Healthcare manometry or reflux monitoring equipment as you work through patient cases, acquire and analyze patient data, and create patient reports.

## Virtual Coaching

Online and in real-time, Diversatek Healthcare Clinical Specialists work with you via screen sharing to provide study-specific data review and report generation coaching for your more difficult studies. Email us at **clinicalsupport@diversatekhc.com** to schedule a one-on-one session.

# Diversatek<sup>™</sup> Healthcare

Advancing GI care by driving science, developing and delivering solutions, and providing unmatched clinical support.

## DiversatekHealthcare.com

P 800.558.6408 or 414.265.7620

### Corporate Headquarters

102 East Keefe Avenue  
Milwaukee, WI 53212 USA

orders@diversatekhc.com  
sales@diversatekhc.com

### Technical Research & Training Center

9150 Commerce Center Circle, Suite 500  
Highlands Ranch, CO 80129 USA

technicalsupport@diversatekhc.com  
clinicalsupport@diversatekhc.com  
clinicaleducation@diversatekhc.com

Learn more from  
your personal product  
representative at  
**800-558-6408**